

COURTNEY THORNE



WIRELESS NURSE CALL SYSTEMS

**DELIVERING CARE TECHNOLOGY
FOR OVER 50 YEARS**



ACCESSORIES



COURTNEY THORNE ACCESSORIES BROCHURE

WHICH COURTNEY THORNE SYSTEM HAVE I GOT?



Altra



Connect 08

All product ID numbers in **blue** or **black** are compatible with this system.




Economy+

All product ID numbers in **purple** or **black** are compatible with this system.


ALTRA & CONNECT CARE RANGE

FULLY COMPLIANT WITH HTM 08-03 GUIDELINES

MAIN DISPLAY PANEL
SERVER ID: 30000



MINI TOUCH DISPLAY
PANEL ID: 30001



CALL BUTTON
ID: 30003



CARDIAC CALL
BUTTON ID:
30005




PEAR PUSH LEAD
ID: 30011



CEILING PULL CORD
ID: 30009



OVER DOOR LIGHT
ID: 30013



ALTRA WEAR ID:
30008



ALTRA GUARD ID:
30010



CORE PRODUCT RANGE

C-T TOUCH

The revolutionary Touch-screen is part of the ultimate Nurse Call System for care homes that require the ability to personalise their call system.

The market leader in resident safety, emergency calls are given priority whilst the optional site map instantly locates the resident, saving vital seconds. Programmable with up to 999 locations in 256 zones, to ensure that the system can be expanded in the future.



ID: **30219**

The C-T Touch gives you the ability to personalise your Nurse Call System, enabling you to safeguard your home by managing your call data for evidencing.

To receive a quotation for adding a C-T Touch to your Courtney Thorne system, please contact our customer service team on 0800 068 7419.

Available to both 08 and Economy+ customers.

DISPLAY PANEL



08 Display Panel

ID: **20007 22521**

The Display Panel features a large, clear, wide viewing angle display with volume control, standby battery and fitted on a swivel bracket.



Economy+ Display Panel

Display Panel Bracket

ID: **22143**



CORE PRODUCT RANGE

ROOM UNITS

The Room Unit allows calls from users and staff. Typical call types include Assistance, Emergency, Visit and Staff Presence.



Room Unit Only
ID: [22835](#) [22823](#)



Room Unit with IP67 Pear Push Lead
(With light switch)
ID: [22825](#) [22831](#) (With light switch) [22821](#)
Typically located in bedrooms, completed with an IP67 rated Pear Push Lead.



Room Unit with Orange Pear Push Lead
ID: [22822](#)
Typically located in bedrooms, completed with an Orange Pear Push Lead.



Room Unit with Pull Cord
ID: [22827](#) [22820](#)
Typically located in bathrooms, completed with Pull Cord.

CORE PRODUCT RANGE

PEAR PUSH LEADS

The Pear Push Leads complement existing Room Units and ensure resident comfort by making it easier for them to call for help if needed.

IP67 Pear Push Lead

ID: 30011 [30012](#) (With light switch)

Designed to withstand complete immersion and is completely dust proof.
Available to both O8 and Economy+ customers. Pear Push with light switch requires Light Switch Relay. Please see ID 22208.



Orange Pear Push Lead
ID: [22117](#)

TRIGGERS

BED & CHAIR WATCH SYSTEMS

The Bed and Chair Watch Systems are designed to assist with residents who may suffer falls and accidents. The systems are simply plugged into the Room unit and will alert staff if the resident gets out of the bed or chair.

BED WATCH SYSTEM

Please note that the solution configuration is system specific for this product. If you are unsure which solution is compatible with your system then please call our customer service team on **0800 068 7419**.

ID: **22251**

Available to both 08 and Economy+ customers.



Timer Unit Connection Lead
ID: 22618

Replacement Power Supply
ID: 22336

Timer Unit
ID: 30350



Wireless Transmitter
ID: 30416

Timer Unit Bracket
ID: 22335



Chair Watch Pad Only
ID: 22182 30411 30412

Available to 08 customers who had their Courtney Thorne system installed post July 2012. Available to Economy+ customers who would like to purchase a spare pad for an existing Chair Watch System.



Chair Watch System
ID: 22256 30414

Available to 08 and Economy+ customers.



Bed Watch Pad Only
ID: 22222 30409

Available to 08 customers who had their Courtney Thorne system installed post July 2012. Available to Economy+ customers who would like to purchase a spare pad for an existing Bed Watch System.



Wireless Bed Watch Pad Only
ID: 30430



Wireless Chair Watch Pad Only
ID: 30431

TRIGGERS

PRESSURE MATS

Pressure Mats are designed to assist with residents who may suffer falls and accidents. The mats simply plug into the Room unit and are placed on the floor. When pressure is applied, an alert is sent to staff.

All pressure mats include a Jack Plug Adaptor, please see page 14 ID:22529.



Dementia Pressure Mat

ID: 30401

The mat is cream in colour instead of the standard black versions; as some people with dementia will not walk on dark spaces on the floor.



Foldable Heavy Duty Pressure Mat

ID: 22790

The mat is fully portable and can be completely folded in half when stored.



Heavy Duty Carpeted Pressure Mat

ID: 22791

This mat has a robust carpet on the top, adding extra fall prevention; it has a heavy non slip backing and is sealed, enabling it to be easily cleaned in case of accidental spillages.



Standard Pressure Mat

ID: 22090
The mat is plastic and is usually recommended for temporary measures, for example respite

WRIST & NECK PENDANTS

08 Wrist and Neck Pendant

ID: 22641

The lightweight pendant comes with a wrist strap or neck anti-ligature lanyard, giving you flexibility in tailoring your approach to care and offering independence to your residents. The pendant is IP67 rated and is both waterproof and dustproof.

Available to 08 customers only.



Economy+ Neck Pendant

ID: 22089

Being easily reprogrammable giving you flexibility in tailoring your approach to care and offering independence to your residents. This lightweight pendant comes with a neck anti-ligature lanyard ensuring wearer safety.

Available to Economy+ customers only.



HEAD SWITCHES

The Head Switches can be used to assist staff with residents who have limited movement or mobility; they can also be the perfect solution for tetraplegics and quadriplegics.

Ping Pong Trigger

ID: 22171

The Ping Pong Trigger is super sensitive to light action and can be activated by subtle movements. The base of the trigger has a universal mounting clip; giving you flexibility with tailoring where to clamp it securely.



TRIGGERS

SPECIALISED TRIGGERS

We offer a number of specialised triggers to further safeguard your care home with evidencing. They are the perfect solution for tailoring your Nurse Call System for residents with individual needs.



Infrared Trigger
ID 30167

The Infrared Trigger has been designed to assist with the problem of residents who are inclined to walk undetected and may suffer falls and accidents. Once installed, the trigger will be activated when the resident breaks the infrared beam.



Staff Attack IR Pendant
ID: **22227**

The Staff Attack Pendant can be carried by staff and can be set off in the event of a dangerous situation arising. Requires charger. Please see page 15 ID 22249 & 22228. Available to 08 customers only.



Enuresis System
ID: **22237**

The Enuresis System sends an alarm to the Nurse Call System upon detection of moisture in the residents bed i.e vomit or urine. Includes a cotton sheet.



Infrared Trigger Bracket
ID: **22221**

The bracket can be used to add extra flexibility with positioning when fixing the Infrared Trigger to the side of the resident's bed.



The Big Button
ID: **22198**

The Big Button can be used by residents with limited dexterity in their hands. It can also assist residents who are experiencing visual deterioration.



Enuresis System Cotton Sheet
ID: **22241**

The Sensing Cotton Sheet has tiny discreet metal fibres threaded into it, which works with the Enuresis Monitor to detect moisture in the residents bed and sends a signal to the Nurse Call System



Epilepsy System
ID: **22237**

Epilepsy Detector comes with Bedside Monitor, Jack Lead, Motion Sensor, Foam, PSU & Rechargeable Battery

PAGING

A Paging System can be integrated with the Nurse Call System to help carers manage multiple tasks at once, by notifying them of a call without the need to visit a touchscreen or display panel. Paging can also be used to reduce ambient sound at night.

If you would like to have a Paging System integrated into your Nurse Call System then please contact our customer service team on **0800 068 7419**.

PAGERS



GEO87ZM Pager
ID: 30139

The pager has a dual line notification screen and uses both characters and digits.

TRANSMITTERS

The Paging Transmitter is part of the Paging System which retransmits the alerts to the Pagers. It is compatible with both the Touchscreen and Display Panel. Please note all new system installations will require an OFCOM licence. If you're unsure please contact us.

X Lite Paging Transmitter
ID: 30174



Pager Rubber Boot for GEO87ZM
ID 22167

A Rubber Boot Holster to protect the GEO85ZM Pager from accidental damage. The rubber boot will include a Screen protector and a Pager strap.



Pager Holster
GEO85ZM
ID: 22266



Pager Battery Cover
GEO85ZM
ID: 22268



GEO14ASM
ID: 22169

Spare Pager Battery Covers for our whole Pager range.

PAGING ACCESSORIES

A selection of spare accessories for our Paging Systems.

CUSTOM PRODUCT RANGE

Light Switch Relay ID: 22208

The light switch relay enables the resident to switch the light on and off from the IP67 pear push lead.



Ceiling Mounted Pull Cord ID: 22191

The Ceiling Mounted Pull Cord is typically located in en-suites / bathrooms and public toilets. Features include an anti-microbial wipe clean cord which is customer replaceable.



Individual Staff Attack IR Triggers Charger ID: 22249



10 x Staff Attack IR Triggers Charger ID: 22228



Re-transmitter / Sounder ID: 22094

The Re-transmitter can be used to boost a weak signal around the home. The product can also be configured into a Sounder which is designed to give an audible indication when a call is activated when a visual Display Panel is not necessary.



Door Contact ID: 22172 22173

The Door Contact is installed onto the door and activates when the door is opened. A key operated timer system gives staff the ability to pass through the door without it alarming.



GEO85ZAAM Pager Charger ID: 22882

Includes a power supply



Over Door Light ID: 22096

The Over Door Light provides an instant visual aid for carers when a resident has alerted for help, it lights up in a specified set of colours to demonstrate the urgency of the call type. Typically installed outside of the room along corridors, it sounds an audio cue and provides a clear visual indication of where the resident is located and when.



Duracell Battery Charger ID: 22232

The charger will charge up to 8x AA & AAA rechargeable batteries and up to 4x C & D rechargeable batteries. It can also charge 1x PP3. Different size cells can be charged at the same time.



Chrome Doorbell ID: 22669

Available to 08 customers only.



Plastic Doorbell ID: 22196

Available to Economy+ customers only.



External Doorbell ID: 22234




REPLACEMENT PARTS

BRACKETS & CLIPS


A selection of Brackets and Clips for all of our Room Units and Pear Push Leads.

Pack of 6 Brackets & 6 Pear Push Lead Clips


ID: 30239




ID: 22705



ID: 22155



ID: 22883



Single Room Unit Bracket

ID: 30019



22386



22185



Pear Push Lead Clip

08 IP67 Pear Push Lead
ID: 22862



Economy+ IP67 Pear
Push Lead
ID: 22577



Economy+ Orange
Pear Push Lead
ID: 22387



ROOM UNIT LABELS

A selection of Room Unit Labels which can be replaced due to gradual wear and tear of the Room Unit buttons, the stickers will help ensure that the buttons are clearly visible to the residents.

08 and Economy+ Labels – Pack of 6

ID: 22276



22866



REPLACEMENT PARTS

ROOM UNIT ACCESSORIES

15m Anti-Microbial Pull Cord ID: 22695

Replacement Pull Cord for the Ceiling Pull Cord. The cord is anti-microbial for extra hygiene making it easy to clean and is waterproof.



Plastic Pull Cord Ring ID: 22874

Replacement Pull Cord Ring for the Pull Cord Room Unit and Ceiling Pull Cord.



Jack Plug Adaptor ID: 22529

The Jack Plug Adaptor can be inserted into one of the auxiliary ports on the Room Units to create additional ports for extra Triggers. It allows carers to tailor the Nurse Call System to use multiple Triggers at the same time.



DOOR CONTACT ACCESSORIES

A selection of spares for the Door Contact.

Door Contact Magnet Only ID: 22146

A replacement Door Contact Magnet.



Door Contact Keys ID: 22148 QTY: Pack of 2

A pack of replacement Door Contact Keys.



AERIALS

Replacement Aerials for the Touchscreen, Panels and Re-transmitters / Sounders.

Whip Aerial

White
ID: 22270



Black
ID: 22696



Dipole Aerial ID: 22264

A more sensitive aerial used to extend signal range.



OTHER ACCESSORIES

LEADS

Jack to Jack Cable ID: 22618

A replacement jack to jack cable which is used on a variety of our Triggers.



NECK PENDANT ACCESSORIES

Neck Pendant Lanyard ID: 22161

Replacement Lanyard for the Economy+ Neck Pendant.



POWER SUPPLIES

A selection of replacement or additional power supplies that cover all the components of the Nurse Call System that need a power supply.

Power Supply for Over Door Lights & Re-transmitter / Sounder ID: 22097



Requires 240v 3pin socket.

STAFF ID TAG ACCESSORIES

Staff ID tag ID: 22231



Staff ID tag – Self Retracting Clip ID: 22244

A spare self-retracting clip for the Staff ID tags.



DOORBELL ACCESSORIES

Door Bell Connection Kit ID: 22265

A Door Bell Connection Kit to assist with integrating a door bell system into the Nurse Call System.



Fallsavers Power Supply Unit ID: 22336

For use with Bed & Chair Watch systems (pages 8 & 9) and Infra Red Trigger (page 11, ID 22209).



OTHER ACCESSORIES

BATTERIES

A selection of replacement batteries for all of our products.

Box of 24 Room
Unit Batteries
ID: 22132



Box of 5 Ceiling Pull Cord / Door
Contact Batteries
ID: 22137



Re-transmitter / Sounder
Battery
ID: 22134



Display Panel
Battery
ID: 22138



Box of 5 Staff ID Fob
Batteries
ID: 22878



Box of 40 Pager
Batteries
ID: 22879



Pager Rechargeable
Batteries
ID: 22207



Economy+ Neck Pendant
Battery
ID: 22135



Box of 5 Wrist and Neck Pendant
Batteries
ID: 22880



TECHNOLOGY FOR THOSE WHO CARE

Courtney Thorne Ltd
Elliott Road
Bournemouth
BH11 8JS
United Kingdom

T: +44(0)1202 573 300

E: info@c-t.co.uk

F: info@c-t.co.uk

W: www.c-t.co.uk

